

Name:

Chart #:

Date:

Snap-On Smile®/Snap It™ Instructions

You have the following treatment(s) done in our office today:

- Snap-On Smile®, Upper arch Snap It™, Upper arch
 Snap-On Smile®, Lower arch Snap It™, Lower arch

Is Snap-On Smile®/Snap It™ like wearing dentures?

The appliance is removable and simply snaps in and snaps out. It is held in place by your natural teeth and does not require adhesives or glue, of any kind.

Insertion & Removal

1. Determine the correct path of insertion; the path of removal is the same as the path of insertion
2. **Always using both hands, start at the end or the back of the appliance.**
3. **Insertion:** apply even pressure until seated completely as demonstrated by your dentist at time of delivery.
4. **Removal:** gently pull the appliance away from the teeth; once the appliance is removed, do not twist, bend, flex, or torque it.

Can I chew gum while wearing Snap-On Smile®/Snap It™?

1. It is not recommended that you chew gum or eat sticky foods.
2. You can eat and drink most food and beverages with it in as you normally would.

Will wearing Snap-On Smile®/Snap It™ damage my natural teeth?

No. Simply follow the use and care instructions.

Will Snap-On Smile®/Snap It™ affect my speech?

1. There will be a temporary adjustment period as you get used to having something new in your mouth. This is very common and most people report that their speech is back to normal within **2 weeks**.
2. Wearing the appliance all day will allow you to return to your normal speech much quicker.

Will Snap-On Smile®/Snap It™ stain?

It does not stain easily, but just like your natural teeth, certain products such as coffee, tea, wine, soy sauce or soda can cause staining.

How to take care of Snap-On Smile®/Snap It™?

1. Clean with a soft bristle toothbrush **without toothpaste** (toothpaste is abrasive).
 2. Clean with **traditional denture soaks**.
 3. Clean with an **ultrasonic cleaner** (purchase online, between \$30 and \$40).
 4. Special cleaning solution & anti-bacterial gel are supplied when you 1st receive appliance, can be reordered through your dentist.
- Instructions: ½ capful of **cleaning solution**: ½ cup of cool water; soak **15 minutes**; water should cover the appliance.

Rinse the appliance thoroughly under running water before using, you may brush it.

Anti-bacterial gel: a small bead into each tooth compartment of the appliance; for best results, apply after each meal.

How long will my Snap-On Smile®/Snap It™ last?

With proper care, your appliance can last from 3 to 5 years.

Limited Warranty

1. It comes with a **12 month** limited warranty against manufacturing defects, from the date the dentist received the appliance.
2. **“One-time Replacement”:** “Denmat makes no warranty with respect to the replacement arch”
3. **Optional Extended Limited Warranty: \$99, two (2) additional years** beyond the standard Limited Warranty, is **valid for 30 days after receipt of appliance.**

DenMat (Manufactory) Customer Service: 1-800-433-6628 www.denmat.com

I am happy with the look, color, size and the bite of my Snap-On Smile®/Snap It™ appliance

Patient's Signature: _____

Legal Guardian's Signature: _____

Relation: _____

Date: _____



Phone: 916-983-9909

whingmingledds@gmail.com

Fax: 916-983-9911